

## TRI-OR™ Surgical Kit & Prosthetics

### Multi Function Orthodontic Implants

#### TRI-OR™ Implant Features

- One-Piece Construction
- Self-Tapping/Drilling
- Made of Grade 5 Titanium
- Cross Slot Bracket Design for 0.025 wire
- 0.7mm Diameter Hole for Round Wire
- 0.7mm Slot for Elastic Band or Connecting Wire

#### TRI-OR™ Implant Sizes

DIAMETER	LENGTH	SKU
1.3mm	6mm	TRI-OR-1.3-6
1.3mm	10mm	TRI-OR-1.3-10
1.8mm	6mm	TRI-OR-1.8-6
1.8mm	10mm	TRI-OR-1.8-10

#### Surgical Instrumentation

HR-TW - Hi-Tec Torque Wrench




HR - Hex Ratchet



ORT - Insertion Tool for TRI-OR Implants



HSD - Hex Straight Driver Handle





**Hi-Tec Implants™** meets and exceeds the highest standards in the field of medical devices: the main approvals, besides many others, are:

**FDA APPROVAL:** Center for Devices and Radiological Health in the US FDA (Food and Drug Administration) Since 1994.

**CE MARK** – After demonstrating compliance with Annex II of Medical Devices Directive 93/42/EEC, entitles us to use CE Marketing on our products.

**ISO 13485:** 2003 – An international standard for quality management of medical devices, Hi-Tec Implants LTD™ meets the requirements of ISO 13485 : 2003 for the design, manufacturing and inspection of dental implants and accessories.

**ISO 9001: 2000** – Certifies that Hi-Tec Implants LTD™ demonstrates compliance of our quality system to meet the requirements of ISO 9001: 2000 (an international standard for quality management system).

**Health Canada Medical Device License and CMDCAS ISO 13485:** 2003 Accredited Since 2005.